LANAFORM®

Massage Armchair

Relax your body

Massage armchair

Fauteuil massant

Massagestoel

MASSAGE ARMCHAIR

EN

User Instructions



Thank you for having chosen the Massage Armchair from LANAFORM®, a high quality massage appliance.

The Massage Armchair is a massage device which provides a deep massage and relaxes your body muscles. In addition, this massage appliance is operated by means of a control panel that is easy to use. LANAFORM® offers you everything you need to profit from the latest advances in the field.

We thank you for your trust.

Names and functions of the components

1. Protective cushion
2. Control panel
3. Arm rest
4. Cushion
5. Support for the legs and foot rest
6. Cushion for the back
7. Air cushion at the arms
8. Jade warmer
9. Mini tray
10. Covering of the back of the seat
11. Opening of the fan
12. Position of the remote control
13. Position of the mini tray
14. Headphone line
15. Headphone
16. MP3 player
17. Connection for the heating jade surface
18. Rear wheels
19. Cables and plug
20. Plug for the headphone
21. Male plug
22. Switch
23. Fuses


Important Safety Warning !

Please read this instruction manual attentively before operating the armchair.

- All use other than shown in this instruction manual is prohibited.
- Please use a supply of electricity suitable for this product.
- Please remove the plug after use or before cleaning the armchair to avoid any injury or damage to the product.
- Do not use this product out of doors.
- Please do not use it if the leather covering has been damaged
- Please do not use it if the ventilating fan of the radiator at the back of the armchair is covered
- Please do not spill anything on this product.
- Please do not use it in conditions of high temperature or in a humid space such as the bathroom, for example.
- Please do not use it in a very dusty room.
- Please do not use it in a room which is not sufficiently ventilated.
- It is preferable to place the massage armchair on a carpet or rug in order to avoid damaging a wooden floor.
- Do not remove or disassemble the MP3 player. Only the USB flash drive can be removed.

Attention: do not move the armchair using the wheels when they are on a wooden floor or an uneven floor in order to avoid damaging the floor. In this case, two persons are required to move the armchair.

This product is not suitable for the following persons:

- Those suffering from osteoporosis.
- Persons suffering from heart ailments or wearing electronic medical equipment such as pacemakers.
- Pregnant women.
- Persons with injuries or suffering from illnesses of the skin.
- Children under twelve or persons who are unable to operate this product.
- Do not use the armchair if you are wet.
- Persons sensitive to heat.

This appliance is not intended for use by persons (including children) whose physical, sensorial or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they are supervised or instructed in the use of the appliance first by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the appliance

Safety

- Make certain that the voltage matches the product specification.
- Do not remove the plug with wet hands.
- Do not let water enter inside this product to avoid electrical shock or causing damage to the product.
- Do not damage the cables or make any modifications to the electrical circuit of this product.
- Use a dry duster to clean this product as well as electrical parts such as the switch and the plug.
- Stay away from this product if it is electrically out of order to avoid any injury should the current return abruptly.
- Stop using this product if you feel ill and go consult a doctor.
- Do not use this product if the outlet is loose.
- Please use it in an intermittent manner. Do not operate this product for hours on end.
- If this product operates without stop for a long period, the thermal protection device will automatically shut it down. The armchair can be used again after a rest of a half hour.

The appliances are not intended for startup by an external timer or a separate remote control system.

If the power cable is damaged, it must be replaced by the manufacturer, an approved repair agent or any qualified person, to avoid any danger.

Maintenance

- This product should be maintained by an approved technician. Users are prohibited from disassembling the armchair themselves.
- Please do not forget to disconnect it after use.
- If you are not going to use this product for a long period of time, fold away the cables and store in a dry place.
- Do not keep this product in an area exposed to high temperature or near an open flame and avoid long periods of exposure to direct sunlight.

Solutions for simple cases of malfunction

- It is quite natural that the motor makes a bit of noise when the armchair is operating.
- If the controls do not work, check to see that the outlet and the plug are well connected and ensure that the switch behind the armchair is lit.
- If the operational time has run out, the product switches off automatically.

Assembling the armchair:

You need to have 2 persons to assemble this armchair.

- 1) Remove the armchair and its accessories from the cartons.
- 2) Place the frame of the armchair on a flat surface
- 3) Place the large mobile protective cushion on the back and secure using the zip (1). We recommend that you avoid moving the cushion behind this mobile cushion to ensure optimal comfort. Then place the small mobile protective cushion on the back and secure using the Velcro strips (2).

Install the armrests of the massage armchair


One person raises the thigh support while the second slides the armrest on the side of the armchair, checking to see that the metal bar slides correctly into the armchair frame. (1)

Attach the screw under the rear of the armrest. (2)





Connect the cables of the air cushion to each case and insert the cases in the extremities of the armrests which you close with protective cover plates

Ensure that the cable behind the protection cushion is connected to the same terminal box as the air cushion.

Connect the headset and the jade heater to the rear of the armchair. (3)


Installation of the control panel and the mini tray

1. Remove the screw that is already installed on the controller stand.


2. Clip the ball that is located in the extremity of the support shelf in the cavity of the controller stand

3


3. Turn the screw in counter clockwise direction

4

4. Remove the 3 screws that are located in the cavity of the armrest

5


5. Insert the metal tip in the armrest

6


6. Screw down the 3 screws in place

7


7. Press down the extremity of the support shelf onto the pre-positioned tip

8


8. Make sure that the 2 holes match

9

9. Screw down the 2 long screws to complete the installation



Assembling the control panel:

Slide the control panel into its support

Assembly of the mini tray:

Place the base of the mini tray in the cavity of the left armrest. Raise the rubber protection and screw down the 3 screws.

Functions and Capabilities

Capabilities

Eliminates fatigue
Encourages blood circulation
Eases muscle fatigue
Relaxes muscle tension
Comforts muscle pains

Functions

This armchair has been designed with a quiet massage mechanism that is activated by four wheels that can rise and descend in the back.

The seat has five air cushions which provide various types of massage using air pressure.


“Magnetic” function: there are 8 magnets (at 300 gauss) both in the area of the buttocks and the area of the foot and calf. The magnets provide a concentrated magnetic power that helps to reduce muscle, joint and rheumatic pain.

Starting the massage

- 1) Press the red button to start the massage armchair (9)
- 2) The mechanism detects the point of massage on your body and automatically adjusts for the position of your shoulders. (If the mechanism for massaging the shoulders is not well adjusted to your position, you can regulate it by using the button «up» and «down.” There are 5 degrees of position for the shoulders.)

 Adjust shoulder (4)	move the mechanism upward to reach the shoulder.
 Adjust shoulder (4)	move the mechanism downward to reach the shoulder.

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The buttons for programming the function A, B, C will light up (9)



Automassage :

It is possible to start a pre-programmed type of massage using buttons 1, 2, 3, 8, 9.

The available pre-programmed massages are: «Comfort,» «Ease,» «Body Care» and «Activate»

The buttons operating the auto massage of the upper body and lower body by air pressure will remain lit. If you press down one of them, it will perform the massage that was entered into the memory previously (1,2,3,8 ou 9)

Name	Button	Massage function	Massage feature	Display status
Characteristic auto massage		Comfort Massage	activate the blood circulation	Auto 01
		Relax Massage	Relax the whole body	Auto 02
		Fitness Massage	Deep massage of the point, relax the ache muscle	Auto 03
		Pressure relieving Massage	Improve the body's fatigue, make the body active	Auto 04
Auto massage		Neck shoulder Massage	Mainly massage the neck, shoulder area	Auto F1
		Back waist Massage	Mainly massage the back waist area	Auto F2
		Lower body pressure massage. Intensity can be adjust	Mainly massage the feet, calf, leg and buttock	Lower body

To change the position of the armchair, use the buttons to adjust the back, legs and angle of recline.

Backrest angle adjustment (button 7)

Backrest raise		Press this button and the backrest will raise slowly, loose the button and it will stop.
Backrest fall		Press this button and the backrest will fall slowly, loose the button and it will stop.

Foot rest angle adjustment (button 7)

Foot rest raise		press this button and the foot rest will raise, loose the button and it will stop.
Foot rest fall		press this button and the foot rest will fall, loose it and it will stop.

Lies angle adjustment (button 7)

Foot rest raise		press this button and the foot rest will raise, loose the button and it will stop.
Foot rest fall		press this button and the foot rest will fall, loose it and it will stop.


Adjust the type of massage

In order to change the position of the armchair, make use of the buttons for adjusting the back, the legs and the angle of recline.

Manual massage of the upper body:

Programmed with five types of massage: Shiatsu, Knocking, Kneading, Flapping, Synchronised Kneading and Flapping.

Manual massage of the upper body (bouton 10)

Function	Button	Description, display status	Massage balls adjustment
Measure	Knocking 		Speed, strength, width
	Shiatsu 		Speed, width, strength
	Kneading 		Speed, strength
	Flapping 		Speed, width, strength
	Kneading & flapping 		Speed, width

You can regulate the speed and width (bouton 11)

In the massage modes kneading, shiatsu and knocking, the space between the two massage balls can be regulated in three width positions: ample, medium and narrow.

Boutton	Display status	Adjust massage balls
Backrest Stretch 	Fixed position (fixed) massage	
	Small area(partial) trip massage	
	Upper body(overall) trip massage	

Speed	SPEED	Speed  activate the blood circulation 5 speed options available
Width	WIDTH	Width  3 width options available

Function: Massage position adjustment (the massage position can be adjusted under fixed and partial modes)

Foot rest raise		Mechanism up
Foot rest fall		Mechanism downt

Manual massage of the lower body (12, 6)

Function	Button	Description, display status		Speciality
Lower body		Lower body air pressure massage	Lowerbody	Two intensity options adjustment
Feet		only foot&calf area air pressure massage	foot&calf	Two intensity options adjustment
buttock		only buttocks air pressure massage	buttocks	Two intensity options adjustment
Intensity		Intensity of the massage by air pressure of the lower body, feet and buttocks.	 intensity	

Other manual functions (15, 13, 5, 14)

Function	Button	Description, display status		Speciality
Vibration		Buttocks vibrating		Two intensity options adjustment in the zones of the buttocks ; strong, weak
Arm		Arm pressure massage foot&calf		Mainly for both arm area
Backrest Heating		Far-infrared heating therapy in the back, you can feel damp heat when you turn it on 3 mins late		The temperature can vary 40°C in 50°C
Feet area Heating		Far-infrared heating for feet area		The maximal temperature can rise to has 40 °C.

The manual function «warm up the jade» (16) Once switched on, the heating of the jade is programmed for 55°. Function of pre-adjustment of the temperature from 40 to 70 degrees, with 5 temperature intervals.


Function	Button	Description, display status		Speciality
Handed heater		Handed heater switch		
		Increase the temperature the highest is 70 degrees		Increase in 5 degrees
		Decrease the temperature the lowest is 40 degrees		Decrease in 5 degrees

MP3 Music function: (17) : you can listen to music whether the armchair is operating or not. The MP3 player has been designed with controls to adjust the volume, listen to the preceding or the next song.

Function	Button	Description, display status	Speciality
Music Play		Music on/off figure 	Music play under boot-strap or awaiting work
		Previous song 	
		Next song 	
		Volume increase 	
		Volume decrease 	

You can enter up to 3 massages in the memory.

Before the end of the massage, press the red button. Frame 8 will flash. Press on one of the letters - A, B or C. The massage is entered in the memory when the words 'memory set on' appear on the screen.


Cut the electricity and stop the massage.

During a massage, press the button 'on/off' two times to stop all massage functions. The back and the footrest will return automatically to their initial position.

After a massage, the armchair will switch off automatically.

Method of downloading MP3 music

- (1) MP3 player
- (2) U-disk front cover
- (3) U-disk

- 1- Disconnect the massage armchair and open the cover of the MP3 player at the back of the armchair.
- 2- Carefully remove the USB stick.
- 3- Insert the USB memory drive into the read/write memory device attached to your computer.
- 4- Download your favourite MP3 formatted songs onto the USB memory drive.
- 5- Remove your USB memory drive from the computer port.
- 6- Insert the USB memory drive into the MP3 player at the back of the armchair and check to see that it functions.
- 7- The MP3 music is updated and you can reconnect the massage armchair to the electricity.

If the MP3 cannot be played, please format the USB memory drive via your computer and start the operation over again.

Product Specifications

Name:	Massage Armchair
Electrical power supply:	110-120V 60Hz 220-240V 50Hz
Maximum power:	300W
Safety structure:	type I
Materials:	PVC, PA, parts made of steel, electrical parts and electronic parts

Limited guarantee

LANAFORM® guarantees that this product is free of all material and manufacturing flaws for a period of two years starting from the date of purchase, except for the points stated below.

The guarantee on this LANAFORM® product does not cover damage caused by improper use or wrong use, accident, the attachment of any unauthorised accessory, any modification made to the product or any circumstance of whatsoever nature that is outside the control of LANAFORM®.

LANAFORM® will not be held responsible for any incidental, consequential or special damages.

All implicit guarantees of the fitness of this product are limited to a period of two years starting from the initial date of purchase.

Upon receipt, LANAFORM® will either repair or replace your appliance and send it back to you. This guarantee is only performed through the LANAFORM® Service Centre. All maintenance activity of this product entrusted to persons other than the LANAFORM® Service Centre cancels the present guarantee.

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SN

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a full page of blank, lined paper. It features approximately 20 evenly spaced horizontal grey lines across its entire width, providing a template for handwriting practice or general note-taking. The margins are consistent on all sides.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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